

Imaging Plate System Portfolio

A film-like workflow that makes going digital easy

Image plate technology provides many practitioners a better compromise between film workflow and digital benefits. Therefore, Carestream Dental offers two imaging plate systems to meet the needs of smallest and largest practices alike—the compact CS 7200 for chairside use and the CS 7600 for multiple users. These systems offer flexible and comfortable plates in the same size as film (without the need for processing); easy image acquisition within seconds; and excellent image quality.

One of the greatest benefits of image plate technology is its similar handling and workflow to analog film. This makes for a smooth transition for practices switching to digital. Plus, an imaging plate system is a great solution to renew digital equipment, which means limited training and less of a learning curve when making the switch.

CS 7200 Imaging Plate System: Your ideal chairside companion

The affordable and ultra-compact CS 7200 digital imaging plate system is an ideal solution for chairside use, offering a cost-effective digital intraoral solution with no compromise to image quality.

- Chairside system designed for daily intraoral indications (size 0, 1, 2)
- Leading image quality with true resolution up to 17lp/mm
- Affordable solution with no compromise to image quality
- Slim, space-saving design
- Fast, streamlined workflow with no clicks required
- Quiet scanning ensures more comfortable exam
- Intuitive workflow prevents mistakes



CS 7600 Imaging Plate System: The intelligent system that's designed for shared users

A digital imaging plate system designed for multiple users, the CS 7600 features a truly automated and secure workflow that delivers outstanding image quality quickly and effortlessly.

- Multi-user system ideal for all intraoral indications (size 0, 1, 2, 3, 4)
- True bitewing and occlusal plate sizes
- Leading image quality with true resolution up to 17lp/mm
- Scan & Go technology provides a truly automated and secure workflow
- Reduces risk of mistakes like mixing up patient images
- Built-in memory prevents loss of data



- Compact unit can be used chairside or shared
- Built-in color monitor provides clear instructions and image previews

Quick Glance at the CS 7200 and CS 7600

	CS 7200	CS 7600
Top line advantages	<ul style="list-style-type: none"> • Best value • Slim, space-saving design • Leading image quality • Quiet scanning 	<ul style="list-style-type: none"> • Scan & Go technology • Automated and secure workflow • Leading image quality • Ideal for all intraoral indications
Use	Chairside / Single user	Shared / Multi-user Chairside / Single user
Max. theoretical resolution	25lp/mm	25lp/mm
Max. true resolution	17lp/mm	17lp/mm
Plate sizes	0-1-2	0-1-2-3-4
Time to first image²	As little as eight seconds	As little as five seconds
FMS workflow	Manual	Automatic
Dimensions (W x H x D)	130 x 270 x 300 mm	237 x 267 x 260 mm

¹ Resolution varies based on the scanning mode and resolution

² Time to first image varies based on the scanning mode and resolution

³ Additional accessories can be ordered separately

To learn more about the CS 7200 or CS 7600, or for more information about any of Carestream Dental's innovative solutions, please call **800.944.6365** or visit **www.carestreamdental.com**.

