

CS 1500



The ideal communication starter

The CS 1500 intraoral camera is the ideal communication tool for any dental practitioner. Delivering precise, true-to-life images with each shot, the CS 1500 camera provides the powerful evidence you need to educate patients and to help them visualise their diagnosis and treatment plan more realistically. Now you can let patients see exactly what you see – so they make smarter decisions, and are more likely to accept treatment.

Industry-leading resolution and focus range

Boasting the industry's highest still image resolution, the CS 1500 camera provides all the functions you need to capture sharp, clear colour images. And with the widest focus range available (1mm to infinity), practitioners can clearly view even the smallest details – making treatment planning and patient education easier and more accurate than ever.

Wireless that works

Available in both wired and wireless configurations, the CS 1500 camera's wireless version provides complete

freedom of movement. Featuring proven Wi-Fi technology, the camera provides faster transmission times, more stability, and less noise than analogue wireless devices – ensuring maximum mobility and reliability.

Effortless image capture

The CS 1500 camera's one-button image capture allows anyone on your team to capture clear, detailed images quickly and easily. Combined with a true autofocus, the camera's unique liquid lens technology works like the human eye to ensure clear, detailed images each and every time.

Flexible, all-purpose design

The CS 1500 camera is a versatile tool for any dental practice. With just one camera, you can capture still images and fluid videos, perform both intra and extraoral exams, and quickly display images on your choice of a PC or video monitor. Meanwhile, USB connectivity allows you to share the docking station between several treatment rooms, which means you can use the camera anytime and anywhere you need it.

Simply integrated

Easily modify, store, and manage your images and videos with the included Carestream Dental's imaging software. Simply press the acquisition button to capture a still image and it's automatically saved to your computer or SD card. Images can then be integrated into patient files, printed, and shared.

Reliable services and support

We provide the guidance and support you need to get the most out of your CS 1500 camera. Our experienced support experts are here to help you integrate your camera into your practice and are available to answer your questions before, during, and after your purchase.



Technical Specifications

Sensor	Micron 1/2.5" CMOS 2592 (H) x 1944 (V)
Video resolution	640 (H) x 480 (V)
Image resolution	1024 (H) x 768 (V)
Focus range	1mm to infinity
Angle of view	90°
Field of view	80°
Focus	Autofocus
Light source	8 White LED array
Video output	USB 2.0; TV-NTSC; TV-PAL; VGA; S-Video
Connection	USB 2.0 high-speed interface
Wireless technology	WiFi - transmission range 10 metres
Workgroup configuration:	Supports up to two groups of ten docking stations

The CS 1500 camera features:

- True autofocus with liquid lens technology
- Best-in-class resolution and focus range
- Reliable Wi-Fi technology
- Auto adjustable white LED lighting for optimised and bright illumination
- One-button simplicity

For more information or to place an order in the UK please call **0800 169 9692**
in Ireland please call **+44 1438 245000**
email **sales.uk.csd@carestream.com**
or visit **www.carestreamdental.co.uk**